



## SUPPLIES NEEDED TO UPLOAD DIABETES DEVICES AT HOME

Device	Software	Necessary items	Automatic upload
ALL glucometers (meters) Accu-Chek Freestyle Lite OneTouch True Metrix	Glooko  <a href="https://my.glooko.com/users/sign_in">https://my.glooko.com/users/sign_in</a>	<ul style="list-style-type: none"> <li>• Glooko account linked to our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ USB cord that connects to your glucometer (please contact our office if you do not have a cord that fits your glucometer)</li> </ul> </li> <li style="text-align: center;"><b>OR</b></li> <li>• Smartphone with Glooko app downloaded on phone               <ul style="list-style-type: none"> <li>○ Special cord that is compatible with both your smartphone and your glucometer or Bluetooth enabled glucometer</li> </ul> </li> </ul>	No
Dexcom G6 CGM (using receiver device)	Clarity  <a href="https://clarity.dexcom.com/">https://clarity.dexcom.com/</a>	<ul style="list-style-type: none"> <li>• Dexcom Clarity account linked to our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ USB cord (charging cable) that connects to your receiver device (please contact our office if you do not have a cord that fits your glucometer)</li> </ul> </li> </ul>	No
Dexcom G6 CGM (using smartphone app)	Dexcom Clarity app	<ul style="list-style-type: none"> <li>• Dexcom Clarity account linked to our office (please contact our office if you are a new user)</li> <li>• Auto uploads - Clarity app downloaded on phone</li> </ul>	Yes
Freestyle Libre CGM (using reader device)	LibreView  <a href="https://www.libreview.com/">https://www.libreview.com/</a>	<ul style="list-style-type: none"> <li>• LibreView account linked to our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ USB cord (charging cable) that connects to your pump (please contact our office if you do not have a cord that fits your pump)</li> </ul> </li> </ul>	No
Freestyle Libre CGM (using smartphone app)	LibreView app	<ul style="list-style-type: none"> <li>• LibreView account linked to our office (please contact our office if you are a new user)</li> <li>• Smartphone with LibreLink app downloaded on phone</li> </ul>	Yes
InPen insulin pen	InPen: Diabetes Management app	<ul style="list-style-type: none"> <li>• Smartphone with InPen: Diabetes Management app downloaded on phone</li> </ul>	No
Medtronic Paradigm and 670G insulin pumps	CareLink  <a href="https://carelink.minimed.com/">https://carelink.minimed.com/</a>	<ul style="list-style-type: none"> <li>• CareLink account linked with our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ CONTOUR NEXT LINK blood glucose meter or CareLink USB</li> </ul> </li> </ul>	No

Device	Software	Necessary items	Automatic upload
Medtronic 770G insulin pump	CareLink <a href="https://carelink.minimed.com/">https://carelink.minimed.com/</a>  <b>OR</b>  MiniMed Mobile app	<ul style="list-style-type: none"> <li>• CareLink account linked with our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ “Blue Adapter” (please contact Medtronic if you do not have adapter)</li> </ul> </li> </ul> <p style="text-align: center;"><b>OR</b></p> <ul style="list-style-type: none"> <li>• Smartphone with MiniMed Mobile app downloaded on phone</li> </ul>	Yes, with phone app
Omnipod Eros insulin pump	Glooko  <a href="https://my.glooko.com/users/sign_in">https://my.glooko.com/users/sign_in</a>	<ul style="list-style-type: none"> <li>• Glooko account linked to our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ USB cord (charging cable) that connects to your pump (please contact our office if you do not have a cord that fits your pump)</li> </ul> </li> </ul> <p style="text-align: center;"><b>OR</b></p> <ul style="list-style-type: none"> <li>• Android Smartphone with Glooko app downloaded on phone Special cord that is compatible with both your smartphone and your PDM (cord can be purchased in-app)</li> </ul>	No
Omnipod DASH insulin pump	Glooko  <a href="https://my.glooko.com/users/sign_in">https://my.glooko.com/users/sign_in</a>  <b>OR</b>  Podder Central app	<ul style="list-style-type: none"> <li>• Glooko account linked to our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ USB cord (charging cable) that connects to your pump (please contact our office if you do not have a cord that fits your pump)</li> </ul> </li> </ul> <p style="text-align: center;"><b>OR</b></p> <ul style="list-style-type: none"> <li>• Android Smartphone with Glooko app downloaded on phone               <ul style="list-style-type: none"> <li>○ Special cord that is compatible with both your smartphone and your PDM (cord can be purchased in-app)</li> </ul> </li> </ul> <p style="text-align: center;"><b>OR</b></p> <ul style="list-style-type: none"> <li>• Auto uploads - smartphone with Glooko app downloaded on phone or computer/laptop with Wi-Fi on</li> </ul>	Yes
Omnipod 5 insulin pump	Glooko <a href="https://my.glooko.com/users/sign_in">https://my.glooko.com/users/sign_in</a>  <b>AND</b>  Podder Central account	<ul style="list-style-type: none"> <li>• Create podder central account when you set-up controller device</li> <li>• Link podder central account to Glooko account (please contact our office if you are a new user to glooko)</li> <li>• Pump will automatically upload data to glooko</li> </ul>	Yes
Tandem Tslim insulin pump	T:connect  <a href="https://tconnect.tandemdiabetes.com/">https://tconnect.tandemdiabetes.com/</a>	<ul style="list-style-type: none"> <li>• T:connect account linked to our office (please contact our office if you are a new user)</li> <li>• Computer/laptop/tablet with USB port               <ul style="list-style-type: none"> <li>○ USB cord (charging cable) that connects to your pump (please contact our office if you do not have a cord that fits your pump)</li> </ul> </li> </ul> <p style="text-align: center;"><b>OR</b></p> <ul style="list-style-type: none"> <li>• Auto uploads - T:connect app downloaded on phone and pump connected to phone via bluetooth</li> </ul>	Yes, with phone app